



Zoom KPM-203D-AH32 / Polishing Machine / Preparing Specimens

Summary

The tester is desgned to obtain specular surfaces of various specimens. It is useful for observing the specular surfaces obtained by a grinder with a microscope.

Features

This machine is designed or manufactured to obtain specular surfaces of te specimens, such as metalic materials.

DETAILED INFORMATION

▶ Specification

• Semi-automatic polisher with semi-auto head mounting

- Plate Size: Φ203

Table 1 plate

- Cleaning method by water spray in bowl
- Motor system for minimizing noise and vibration
- Rotation method by plate rotating toe
- Central pressurization system by air cylinder

• Power : AC220V, 60Hz

• Size : 750(W) * 650(D) * 900(H) mm

▶ Application

Molded specimen, products that require gloss