

Distinctive "performance" and "design"
The leader of 4th-generation digital testers

KMT



► Summary

It is to form specimens to facilitate the hardness and polishing works of cut of various materials.

► Features

It is to form specimens to facilitate the hardness and polishing works of cut of metallic and various materials.



Zoom

KMP-130-2 / Mounting Press / Preparing Specimens

▣ DETAILED INFORMATION

► Specification

- Power : AC 220V, 50 / 60 Hz
- Mold Production Time : 6 minutes of shortest heating
7 minutes of cooling
- Mold Size : Basic 32, 40mm (Selectable)
Customizable
- Transfer Distance : 130 mm
- Easy-to-open quick lever system
- Temperature, Maintenance, Cooling are controlled by program

- Equipped with internal anti-dust and anti-oil cover
- Cylinder to be descended and melody warning when job is completed
- Size : 405(W) * 500(D) * 430(H) mm

► **Application**

Easy to make specimens and polish