

# **ELECTRONIC SPIROMETER TRAINER Model Number : GOTT-ESMT-01**





### **DESCRIPTION**

GOTT-ESMT-01 is designed to demonstrate the operation, theory, repair and servicing of Electronic Spirometer equipment that used at clinics & hospitals. The unit is given the fault diagnosis so as to facilitate the teaching of trouble shooting techniques.

A fault simulator is provided with circuit block diagram of the device for the students to study theoretically and experimentally in order to understand the operation of the device.

The fault is activated by button that disables or create malfunction of certain circuit or components so that the device become faulty.

Based on the symptoms and the signal at the test points, the students then do the troubleshooting exercise to find the faulty circuit or components. After locating the failure, the fault can be reset, and then the device work normally.

#### **FEATURES**

- Real Electronic Spirometer unit
- · High quality industrially accepted component
- Bold color identification of circuits
- Designed to simulate common faults
- Detailed instruction & experiment manual

#### **SPECIFICATIONS**

- Flow rate-volume chart, volume-time chart display
- Data memory, delete, upload and review
- Trend chart display
- Automatic power off when there is no operation in one minute
- Circuit Block Diagram of Electronic Spirometer
- Electrical Test Points
- Fault Simulation

## **EXPERIMENT TOPICS**

- Understanding Electronic Spirometer Unit
- Block Diagram Electronic Spirometer Explanation
- Operating The Electronic Spirometer Training Unit
- Troubleshooting of Rectifier

- Troubleshooting of Microcontroller
- Troubleshooting of Control Panel
- Troubleshooting of Alarm
- Troubleshooting of Signal Conditioning

# Manuals:

- (1) All manuals are written in English
- (2) Model Answer
- (3) Teaching Manuals

## **General Terms:**

- (1) Accessories will be provided where applicable.
- (2) Manual & Training will be provided where applicable.
- (3) Design & specifications are subject to change without notice.
- (4) We reserve the right to discontinue the manufacturing of any product.

# Warranty:

2 years

#### **ORDERING INFORMATION:**

ITEM	MODEL NUMBER	CODE
ELECTRONIC SPIROMETER TRAINER	GOTT-ESMT-01	200-408